CHAGAS DISEASE

SD BIOLINE Chagas Ab Rapid TEST

Simple - Easy - Accurate - Rapid

- ICA rapid test for the detection of antibodies to *Trypanosoma cruzi*
- Effective test for acute & chronic Chagas infection
- Specimen : Whole blood or Serum/Plasma (100µl)
- Fast test result : 15 minutes

Chagas disease, or American trypanosomiasis, is caused by the parasite *Trypanosoma cruzi*. Infection is most commonly acquired through contact with the feces of an infected triatomine bug (or "kissing bug"), a blood-sucking insect that feeds on humans and animals.

■ The fact about Chagas disease

- Affect 16–18 million people as of 2008.
- Some 100 million at risk of acquiring the disease. (25% of the Latin American population)
- Kill around 50,000 people annually.
- Present in 18 countries on the American continent. (from the southern United States to southern Argentina).

■ Why is it important to detect antibody of Chagas disease?

The diagnosis of Chagas disease can be made by observation of the parasite in a blood smear by microscopic examination. However, a blood smear works well only in the acute phase of infection when parasites are seen circulating in blood. Diagnosis of chronic Chagas disease is made after consideration of the patient's clinical findings, as well as by the likelihood of being infected, such as having lived in an endemic country. **Therefore, to diagnose chronic Chagas, the serological test for antibody to the parasite is preferable instead of direct parasite detection by microscopic examination.**

And Diagnosis is generally made by testing with at least two different serologic tests.

■ Performance SD BIOLINE Chagas Ab Rapid

		Commercial ELISA		Total
		Positive	Negative	
SD Chagas Ab Rapid	Positive	139 *	0	139
	Negative	1	140	141
Total		140 *	140	280
SD Sensitivity		99.2% (139/140)		
SD Specificity		100% (140/140)		

Sensitivity : 99.2%(139/140)

Specificity : 100%(140/140)

* Positive samples were selected from National Reference Laboratory for Chagas and Leishmaniasis in Honduras.

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