# SD BIO LINE Cholera Ag O1/O139

## Explanation of the test

Cholera is an acute intestinal disease with watery diarrhea, vomiting, high dehydration, acidosis, and circulation disorders. The causative agent of cholera is a gram-negative bacterium, *Vibrio cholerae* and this is generally transmitted to humans via contaminated water and food. Because of the clinical and epidemiological severity of cholera, it is critical to determine as quickly as possible the presence of *V. cholerae* O1/O139 in clinical specimens, water and food so that appropriate monitoring and preventive measures can be undertaken by public health authorities.

Conventional culture methods currently available for detection of *V. cholerae* O1/ O139 are time-consuming and expensive and lack the necessary sensitivity, particularly for clinical samples from patients previously treated with antibiotics and/or traditional home remedies or for poorly handled samples. Polyclonal antibody-based coagglutination test for *V. cholerae* O1/ O139 requires cross-absorption with staphylococci. The PCR has been reported to be a reliable and sensitive method for the detection of *V. cholerae* O1/ O139. However, PCR requires trained personnel and sophisticated equipment, which are not generally available in all areas of the world, particularly where cholera is a serious and recurrent problem.

The SD BIOLINE Cholera Ag O1/O139 kit is a rapid, qualitative test for the detection of *V. cholerae* O1/O139 in human fecal specimens. The SD BIOLINE Cholera Ag O1/O139 test contains a membrane strip, which is pre-coated with mouse monoclonal anti-*V. cholerae* O1 antibody on test line 1 (1) region and with mouse monoclonal anti-*V. cholerae* O139 antibody on test line 2 (2) region on the surface of the device. These lines in result window are not visible before applying any samples. The Control Line is used for procedural control. Control Line should always appear if the test procedure is performed properly and the test reagents of control line are working. A purple Test Line will be visible in the result window if *V. cholerae* O1/ O139 antigen is present in the sample.

When a specimen is added to the test, *V. cholerae* O1/O139 antigen in the specimen reacts with colloidal-gold-labeled *V. cholerae* O1/ O139 specific antibodies and forms a complex of antigen-antibody colloidal gold conjugates. As this complex migrates along the length of the result window by capillary action, the complex is captured by the *V. cholerae* O1/ O139 specific monoclonal antibodies immobilized in the test line across the result window and generates a colored line. In the absence of *V. cholerae* O1/ O139 in specimens, a complex is not formed and no colored line appears in the result window of test device.

## Materials provided

1. Test device individually foil pouched with a desiccant
2. Sample collection tube
3. Assay diluent (25 ml/ vial)
4. Sample collection swab
5. Disposable dropper
6. Instruction for use

## Kit storage and stability

1. The SD BIOLINE Cholera Ag O1/O139 test kit should be stored at 1~30°C.
2. The test device is sensitive to humidity as well as to heat.
3. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use it beyond the expiration.
5. Do not store at refrigerator. Do not freeze.
6. The shelf-life of the kit is as indicated on the outer package.

## Specimen collection, Preparation, Storage and Precaution

1. Sample collection and Preparation
  - To take a portion of feces (about 50mg), insert the sterile swab into a stool sample that presents the most secretion under visual inspection.
  - Insert the swab into the sample collection tube containing assay diluent.
  - Swirl the swab at least 10 times until the sample has been dissolved into the assay diluent and discard the swab while squeezing the swab against the wall of tube, replace the cap.
  - Specimen extracted in assay diluent may be stored at 2 ~ 8°C for up to 1 week prior to testing.
2. Specimen transport and storage
  - Specimen should be tested as soon as possible after collection. Do not use any kind of transport media to store or transport specimens
  - Faecal sample may be stored refrigerated (2 ~ 8°C) for 72 hours. If longer storage is required, freezing at -20°C is recommended.
3. Precaution
  - Faecal specimens for direct test should be collected into containers that do not contain media, preservatives, as all of these additives may interfere with SD BIOLINE Cholera Ag test.

## Warnings & Precautions

1. Precautions
  - SD BIOLINE Cholera Ag O1/O139 test should be stored at room temperature.
  - The test device is sensitive to humidity as well as to heat.
  - Perform the test immediately after removing the test device from the foil pouch.
  - Do not use the test kit if the pouch is damaged or the seal is broken.
2. Warnings
  - For in vitro diagnostic use only. Do Not Re-use test device.
  - The instruction must be followed exactly to get accurate results. Anyone performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
  - Do not eat or smoke while handling specimens.
  - Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
  - Avoid splashing or aerosol formation.
  - Clean up spills thoroughly using an appropriate disinfectant.
  - Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.
  - Do not mix and interchange different specimen.
  - Care should be taken to avoid contamination of the end of sample collection tube when dropping assay diluent into sample well.
  - SD BIOLINE Cholera Ag O1/O139 assay diluent contains a proprietary anti-microbial agent which presents no hazard to the user if normal laboratory safety precautions are followed.

## Procedure of the Test (Refer to figure)

[Extraction Procedure : preparation of extracted sample]

1. Allow test device and extracted sample to room temperature prior to testing.
2. Take assay diluent up to the Fill Line as shown in Figure. And then, transfer into the sample collection tube.
3. Repeat 2)
4. Take a portion of feces(about 50mg) from a stool sample with the sample collection swab.
5. Insert the swab into the sample collection tube containing assay diluent.

6. Swirl the swab at least 10 times until the sample has been dissolved into the assay diluent and discard the swab while squeezing the swab against the wall of tube.
  7. Assemble dropping cap on the sample collection tube.
- [Test Procedure]
1. Remove the test device from the foil pouch, and place it on a flat, dry surface.
  2. Add 3~4 drops (about 80~100µl) into the sample well (S) of the test device.
  3. As the test begins to work, you will see purple color move across the result window in the center of the test device.
  4. Interpret test results at 10~20 minutes. Do not interpret test result after 20 minutes.

## Interpretation of the test (Refer to figure)

1. A color band will appear in the left section of the result window to show that the test is working properly. This band is control line (C).
2. Color bands will appear in the middle and right section of the result window. These bands are test line 2 and test line 1 (2, 1).

**Negative Result :** The presence of only control band (C) within the result window indicates a negative result.

**Positive Result :**

1. The presence of two lines as control line (C) and test line 1 (1) within the result window indicates a positive result for *V. cholerae* O1 antigen.
2. The presence of two lines as control line (C) and test line 2 (2) within the result window indicates a positive result for *V. cholerae* O139 antigen.
3. The presence of three lines as control line (C), test line 1 (1) and test line 2 (2) within the result window indicates a positive result for *V. cholerae* O1 and O139 antigen.

**Invalid Result :** If the control band (C) is not visible within the result window after performing the test, the result is considered invalid. Some causes of invalid results are because of not following the directions correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be re-tested using a new test kit.

## Limitation of the test

1. A negative result does not exclude the possibility of *V. cholerae* O1 and/or O139 infection in a patient. Failure to detect *V. cholerae* O1 and/or O139 may be a result of factors such as collection of specimen at an improper time in the disease when few bacteria are present and improper sampling or handling of the specimen.
2. The SD BIOLINE Cholera Ag O1/ O139 test kit detects *V. cholerae* O1 and/or O139 in human stool specimens.
3. A positive result does not preclude the presence of other enteric pathogens. While the relationship between cholera and gastroenteritis is well established, concurrent infection with other microbial pathogens is possible. Additional microbiological tests should be performed in parallel with SD BIOLINE Cholera Ag O1/ O139 test kit in order to exclude other possible causes of the illness.
4. Test results should be interpreted in conjunction with information available from epidemiological studies, clinical assessment of the patient and other diagnostic procedures.

## Internal Quality Control

The SD BIOLINE Cholera Ag O1/ O139 test device has "Test Line 1", "Test Line 2", and "Control Line" on the surface of the cassette. All the Test Lines and Control Line in result window are not visible before applying any samples. The Control Line is used for procedural control. Control line should always appear if the test procedure is performed properly and the test reagents of control line are working.

## Performance Characteristics

The performance of SD BIOLINE Cholera Ag O1/O139 was confirmed by culture and PCR method. The results are summarized in the following tables.

### 1. Cholera O1 Ag Sensitivity and Specificity

		SD BIOLINE Cholera Ag O1/O139 Rapid		Commercial Cholera Ag O1/O139 Rapid test	
		Cholera O1 Ag Pos.	Cholera O1 Ag Neg.	Cholera O1 Ag Pos.	Cholera O1 Ag Neg.
		<b>Confirmed by Culture and PCR method</b>	<i>V. Cholera</i> O1 Ag Pos. 152 specimens	145	7
	<i>V. Cholera</i> O1 Ag Neg. 320 specimens	19	301	33	287
	Sensitivity	95.4% (145/152)		94.1% (143/152)	
	Specificity	94.1% (301/320)		89.7% (287/320)	

### 2. Cholera O139 Ag Sensitivity and Specificity

		SD BIOLINE Cholera Ag O1/O139 Rapid		Commercial Cholera Ag O1/O139 Rapid test	
		Cholera O1 Ag Pos.	Cholera O1 Ag Neg.	Cholera O1 Ag Pos.	Cholera O1 Ag Neg.
		<b>Confirmed by Culture and PCR method</b>	<i>V. Cholera</i> O139 Ag Pos. 120 specimens	119	1
	<i>V. Cholera</i> O139 Ag Neg. 352 specimens	7	345	15	337
	Sensitivity	99.2% (119/120)		98.3% (118/120)	
	Specificity	98.0% (345/352)		95.7% (337/352)	

## Bibliography of suggested reading

1. Colwell, R.R et al. Development and evaluation of a rapid, simple, sensitive, monoclonal antibody-based coagglutination test for direct detection of *Vibrio cholerae* O1. 1992. FEMS Microbiol. Lett. 97:215-220.
2. Fields, P.I et al. Use of polymerase chain reaction for detection of toxigenic *Vibrio cholerae* O1 strains from the latin American cholera epidemic. 1992. J. Clin. Microbiol. 30:2118-2121.
3. Detection of *Vibrio cholerae* O1 in the aquatic environment by fluorescent-monoclonal antibody and culture methods. 1990. Appl. Environ. Microbiol. 56:2370-2373
4. Environmental *Vibrio cholerae* O139 may be the progenitor of outbreak of cholera in coastal area of Orissa, eastern India, 2000: molecular evidence. Khuntia HK, Pal BB, Meher PK, Chhotray GP. Am J Trop Med Hyg. 2008 May;78(5):819-22. No abstract available.
5. Detection and identification of enterohemorrhagic *Escherichia coli* O157:H7 and *Vibrio cholerae* O139 using oligonucleotide microarray. Jin DZ, Xu XJ, Chen SH, Wen SY, Ma XE, Zhang Z, Lin F, Wang SQ. Infect Agent Cancer. 2007 Dec 23;2:23.

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- Colwell, R.R et al. Development and evaluation of a rapid, simple, sensitive, monoclonal antibody-based coagglutination test for direct detection of *Vibrio cholerae* 01, 1992, FEMS Microbiol Lett. 97:215-220.
- Fields, P. I et al. Use of polymerase chain reaction for detection of toxigenic *Vibrio cholera* 01 strains from the Latin American cholera epidemic, 1992. J. Clin. Microbiol 30: 2118-2121.
- Detection of *Vibrio cholera* 01 in the aquatic environment by fluorescent-monooclonal antibody and culture methods. 1990. Appl. Environ. Microbiol. 56: 2370-2373.
- Environmental *Vibrio cholera* 0139 may be the progenitor of outbreak of cholera in coastal areas of Orissa, eastern India, 2000: molecular evidence, Khuntia HK, Pal BB, Meher PK, Chhotray GP. Am J Trop Med Hyg. 2008 May; 78(5):819-22. No abstract available.
- Detection and identification of enterohemorrhagic *Escherichia coli* O157:H7 and *Vibrio cholera* 01439 using oligonucleotide microarray. Jin Dz, Xu XJ, Chen SH, Wen SY, Ma Xe, Zhang Z, Lin F, Wang SQ. Infect Agent Cancer. 2007 Dec 23; 2:23

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